

imp De Mayo V² E™ Knee Positioner

Versatility with V²E™

De Mayo V² E™ Knee Positioner



◆ When utilizing the Sterile Extension, the surgeon can stand between the legs of the patient for greater access approaching the medial compartment.



De Mayo V²™ Knee Positioner




Illustrated with the De Mayo Universal Distractor®

- ◆ Easier access to the medial compartment improves surgeon posture.
- ◆ Surgeon has an option of either the Extension “off the table” or the Standard “on the table” capabilities.
- ◆ New Carriage with improved holding power.



◆ Varus tilt control allows varying degrees of lateral movement by pushing the slide bar in increments from open to close.





◆ Removable “transformable” carriage facilitates cleaning.



Sterile Gel-infused
memory foam *imp*
Boot Pad &
Cohesive Wrap

Cat. No. 803-GP-10

Ordering Information

Description	Cat. No.
De Mayo V2 E™ Knee Positioner® System	114-SYS KIT
INCLUDES: De Mayo V2™/V2 E™ Knee Positioner® 30" Baseplate w/ Removable Carriage & Single Bar Lock Base	114-BP
De Mayo V2™/V2 E™ Knee Positioner® Clamp	114-CL
De Mayo V2 E™ Knee Positioner® Sterile Extension	114-EXT
Aluminum Boot w/Distractor Block	803-ABD
De Mayo V2 E™ Knee Positioner® System	114-S-SYS KIT
INCLUDES: De Mayo V2™/V2 E™ Knee Positioner® 25" Baseplate w/ Removable Carriage & Single Bar Lock Base	114-S BP
De Mayo V2™/V2 E™ Knee Positioner® Clamp	114-CL
De Mayo V2 E™ Knee Positioner® Sterile Extension - Short	114-S EXT
Aluminum Boot w/Distractor Block	803-ABD
De Mayo V2 E™ Knee Positioner® System	114-907-SYS KIT
INCLUDES: De Mayo V2™/V2 E™ Knee Positioner® 30" Baseplate w/ Removable Carriage & Single Bar Lock Base	114-BP
De Mayo V2™/V2 E™ Knee Positioner® Clamp	114-CL
De Mayo V2 E™ Knee Positioner® Sterile Extension	114-EXT
Aluminum Boot w/Distractor Block	803-ABD
De Mayo Universal Distractor®	907
De Mayo V2 E™ Knee Positioner® System	114-S-907-SYS KIT
INCLUDES: De Mayo V2™/V2 E™ Knee Positioner® 25" Baseplate w/ Removable Carriage & Single Bar Lock Base	114-S BP
De Mayo V2™/V2 E™ Knee Positioner® Clamp	114-CL
De Mayo V2 E™ Knee Positioner® Sterile Extension - Short	114-S EXT
Aluminum Boot w/Distractor Block	803-ABD
De Mayo Universal Distractor®	907



De Mayo V2 E™ Knee Positioner is a trademark of Innovative Medical Products, Inc.
U.S. Patent No. 8,132,278

V2™/V2 E™ Clamp is a trademark of Innovative Medical Products, Inc.
U.S. Patent No. 7,003,827

De Mayo Universal Distractor® U.S. Patent No. 7,380,299

Proprietary Gel Foam Boot Pad property of Innovative Medical Products, Inc.

© October 2017 Innovative Medical Products, Inc.

ALL RIGHTS RESERVED

imp® Innovative
Medical Products, Inc.

87 Spring Lane, PO Box 8028
Plainville, CT 06062

PH: (860) 793-0391 FAX: (866) 459-1805

Toll Free PH: 1-800-467-4944

www.impmedical.com

