

De Mayo Universal Distractor®

New Technique for Femoral Distraction

Improves vision and access to posterior compartment

Clear field of vision allows easier, unobstructed access

Reduce procedure time for Unicompartamental, MIS Techniques, Arthroscopies, ACL and TKA surgeries

Pressure Protector Pad® protects the popliteal area

Requires the De Mayo Knee Positioner® that features a single locking lever for precise control of flexion/extension, tilt & rotation

Broad distribution of pressure with sterile protector pad that will not rotate or roll

U.S. Patent No. 8,048,082

Tension release lever to fine-tune spacing or disengage distraction

U.S. Patent No. 7,380,299

U.S. Patent No. 7,003,827

Hand control for finite adjustment of distraction

PATENT APPLIED FOR

A Solution for Tight Spaces

The De Mayo Universal Distractor® applies finite pressure to the femur outside the knee joint. It allows an unobstructed view into the posterior capsule without the use of lamina spreaders, or a second pair of hands to apply traction. Maximum view allows for easy insertion of instruments, trials, and implants while performing Uni, TKA, ACL and Tibial Plateau surgeries.



De Mayo Universal Distractor®



Surgical Photo courtesy of Robert Kennon, M.D., Waterbury Hospital, Waterbury CT

◇ Improves access to the posterior aspect of the joint for removal of osteophytes, resected bone, bone cement, and cautery applications

◇ Improves MIS techniques for component insertion and prevents component damage

◇ Opens medial and lateral compartments without the use of a “stress” post

◇ Surgeon controlled distraction of the compartments reducing potential over stress to the patient’s soft tissue



Surgical Photo courtesy of John Keggi, M.D., Waterbury Hospital, Waterbury CT

Ordering Information

Description

Description	Cat. No.
De Mayo Universal Distractor® (requires Boot w/Distractor Block)	907
Pressure Protector Pad Set (incl. 294-P Boot Pad, Distractor Pad & Cohesive Bandage)	907-P
De Mayo Knee Positioner® with Universal Distractor® 30”	803-907-ABD
De Mayo Knee Positioner® with Universal Distractor® 25”	803-S-907-ABD

Contact IMP for additional ordering options

imp® *Innovative*
Medical Products, Inc.

A GLOBAL LEADER IN PATIENT POSITIONING

87 Spring Lane, Plainville, CT 06062

PH: 860-793-0391 FAX: 866-459-1805

Toll Free PH: 800-467-4944

www.impmedical.com

De Mayo Universal Distractor® is a registered trademark of Innovative Medical Products, Inc.
De Mayo Universal Distractor® - U.S. Patent No. 8,048,082
Pressure Protector Pad® is a registered trademark of Innovative Medical Products, Inc. -PATENT APPLIED FOR
De Mayo Knee Positioner® is a registered trademark of Innovative Medical Products, Inc.
De Mayo Knee Positioner® - U.S. Patent No. 7,380,299
Single Lever Clamp - U.S. Patent No. 7,003,827
© 2013 Innovative Medical Products, Inc.
All rights reserved 01/2013

CAT. No. 907-P
INCLUDES:

STERILE

Ⓜ SINGLE USE



LATEX FREE FOAM