



**De Mayo Universal Distractor® provides surgeons two important benefits:
precise distraction of the knee during surgery,
and an unobstructed view of the surgical site**

PLAINVILLE, Conn. (August 6, 2018) - Designed by Innovative Medical Products, the De Mayo Universal Distractor® not only enables surgeons to deliver precise, finite distraction of the knee during surgery but also lets surgeons see into the posterior regions of the joint, providing a completely unobstructed view of the operative site.

Used in conjunction with the De Mayo Knee Positioner®, the De Mayo Universal Distractor® is a sterile, external device that only works with the IMP boot and attaches to the back of the boot. It allows the surgeon to independently control the distracting of the knee joint by applying finite pressure to separate the joint space by means of a handle on the device that incrementally pumps the required tension.

During Arthroscopic procedures, the IMP solution eliminates the need for the surgeon to achieve distraction by applying pressure to the joint with the aid of a stress post, which sometimes creates post-operative pain.

Reducing procedure time for Unicompartmental and TKA surgeries, the IMP distractor also eliminates the need for lamina spreaders, the use of a bone hook, or having an assistant manually pull on the femur at the thigh throughout the surgical procedure. The IMP solution is like having an extra pair of hands in the OR.

Pressure Protector Pad®

The De Mayo Universal Distractor® comes with a sterile, single-use, latex-free Pressure Protector Pad® set that protects the popliteal area of the knee from direct pressure to the popliteal nerve and soft tissue. Designed with an inner core of dense, non-collapsible plastic material covering the square tube portion of the distractor, the pad cannot roll or rotate, providing a wide distribution of forces protecting the delicate structures of the knee.

Customers:

For more information about the De Mayo Universal Distractor®, please contact Innovative Medical Products, 87 Spring Lane, P.O. Box 8028, Plainville, CT 06062; phone: 860-793-0391 or toll free: 800-467-4944; or visit www.impmedical.com.

About Innovative Medical Products

Innovative Medical Products Inc.® designs, manufactures, and brings to market patient positioning devices for healthcare and sports related surgeries. With more than 40 years of experience in supporting hospitals,

surgery centers, and clinics, IMP provides positioning products that stabilize and protect patients during navigated, computer assisted and robotic surgeries. IMP devices save hospital and surgery center costs by minimizing staff, decreasing liability and improving the overall patient experience. IMP manufactures, patents and trademarks products to Federal and International standards for quality and safety for both the patient and the healthcare institution. For more information, visit: www.impmedical.com

#

Media contact: Harry McBrien, Maier True Communication; 860-470-3127; harry@maier.com