

# Reznik Universal Shoulder Positioner™

Designed for shoulder arthroscopies in the lateral position.

3 planes of adjustment~ Rotation, Horizontal, Vertical

#1 - 360° rotational adjustment

#2 - Horizontal plane with stop @ 90°

#3 - Vertical plane with a stop to prevent components from falling onto patient



Provides controlled distraction of the shoulder.

*Phase 4 Gel™ Splint* increases comfort, protects skin, controls internal & external rotation while holding the arm securely in place.

“Whether it is position, rotation control or flexibility every surgeon’s set up is possible with this simplistic universal light weight design.”

The Reznik Universal Shoulder Positioner™ can be utilized in either the lateral position, Beach Chair and for elbow arthroscopies.

“Solves a wide array of surgical needs with one device. The system represents an economical solution to the challenge of multi-surgeon, multi-surgical technique traction requirements.”



Surgical photo courtesy of Alan Reznik, M.D. New Haven, CT



## Order Information

Description

Cat. #:

Reznik Universal Shoulder Positioner™ System  
[includes: Reznik Universal Positioner™, imp® rail clamp]

908 KIT

Phase 4 Gel™ Splint ~ STERILE  
[10/ case]

509-P4

Phase 4 Gel™ Splint ~ NON STERILE  
[10/ case]

509-P4-NSP



# Tri-Pull™ Secure Shoulder Solution

A complete system for greater access

Tri-Pull™ Secure Shoulder Solution combines three proven systems to secure and distract the shoulder precisely.

Featuring the patented **Reznik Universal Shoulder Positioner™** with three planes of adjustment.



**Phase 4 Gel™ Splint** adheres to arm for secure positioning and features a natural wrist contour for increased comfort and patient safety.



**De Mayo RoTractor®** delivers complete patient control in the sterile field, holding fixed positions of rotation and allowing distraction of the shoulder or elbow.

Ideal for shoulder or elbow arthroscopies~

## Order Information

Description	Cat. #:
Tri-Pull™ Secure Shoulder Solution [ includes: Reznik Universal Shoulder Positioner™, De Mayo RoTractor®, imp® rail clamp]	908 SYS KIT
Phase 4 Gel™ Gel Splint ~ STERILE [10/ case]	509-P4
Phase 4 Gel™ Splint ~ NON STERILE [10/ case]	509-P4-NSP



Reznik Universal Shoulder Positioner™ is a trademark of Innovative Medical Products, Inc.  
Reznik Universal Shoulder Positioner™ U.S. Patent No. 7,569,024  
De Mayo RoTractor® is a registered trademark of Innovative Medical Products, Inc.  
De Mayo RoTractor® - PATENT APPLIED FOR  
imp® Rail Clamp is a trademark of Innovative Medical Products, Inc.  
Phase 4 Gel™ Splint is a trademark of Innovative Medical Products, Inc.  
Phase 4 Gel™ Splint- PATENT APPLIED FOR  
© Sept 2013 Innovative Medical Products, Inc.  
ALL RIGHTS RESERVED

 **ism** Innovative Sports Medicine™  
A DIVISION OF INNOVATIVE MEDICAL PRODUCTS, INC.  
87 Spring Lane, PO Box 8028  
Plainville, CT 06062  
PH: (860) 793-0391 FAX: (866) 459-1805  
Toll Free PH: 1-800-467-4944  
www.impmedical.com