



# De Mayo RoTractor®



Adapts to most pole style positioners.



De Mayo RoTractor® holds patient's arm in rotation while maintaining distraction of the shoulder.



○ Rotation can be "locked" by tightening locking knob.

} RoTractor® can be autoclaved for sterile use.



Rotation is achieved by simply rotating the patient's arm and is maintained by the RoTractor's uniquely designed gear.

RoTractor® has 3 connector posts matching the hole configuration of the Phase 4 Gel™ Splint for maximum security.



Phase 4 Gel™ Splint



## Ordering Information

<u>Description</u>	<u>Cat. No.</u>
De Mayo RoTractor®	# 011102
De Mayo RoTractor® Adapter - Arthrex	# 011104
De Mayo RoTractor® Adapter - Reznik	# 021204
Phase 4 Gel™ Splint~ STERILE [10/case]	# 509-P4
Phase 4 Gel™ Splint ~ NON STERILE [10/case]	# 509-P4-NSP

**ism** Innovative Sports Medicine®

A DIVISION OF INNOVATIVE MEDICAL PRODUCTS, INC.

87 Spring Lane, PO Box 8028

Plainville, CT 06062

PH: (860) 793-0391 FAX: (866) 459-1805

Toll Free PH: 1-800-467-4944

www.impmedical.com

