

SuperBump™

“ A Bump that holds its own ”

When simplicity is all you need for a knee procedure.



Photo courtesy of Steven B. Haas, M.D., NY



Positioned on top of the drapes at any location for knee flexion.



Contour of the SuperBump™ cradles the leg during closure.

IMP...THE OPERATIVE WORD
IN PATIENT POSITIONING

SuperBump™



STERILE

STERILE EO


SuperBump™ is packaged sterile and accepts two five pound sterilizable stainless steel weights.



② LATEX FREE FOAM

Eliminates cross contamination by guaranteed sterility.

Weights come with a sterilization case for autoclaving and storage.



Ordering Information

<i>Cat. No.</i>	<i>Description</i>
0911 SYS KIT	includes: 2- 091102 SuperBump™ Weights 1- 091102-CASE Sterilization Case
091101	SuperBump™ 5/case



imp® **Innovative**
Medical Products, Inc.

87 Spring Lane, PO Box 8028
Plainville, CT 06062

PH: (860) 793-0391 FAX: (866) 459-1805
Toll Free PH: 1-800-467-4944
www.impmedical.com