

Arm Straps allows access for leads and IV's by Anesthesiology.
Prevents potential neurological impairment caused by sheet tucking.
Designed to meet AORN recommendations.

Estape TrenMAX[®]

Trendelenburg Pad

The Sticky Pad™

Providing superior holding power and clinical use since 2009




Single-use Pad prevents potential cross contamination

TrenMAX[®] prevents nerve damage by eliminating pressure on the brachial plexus

TrenMAX[®] doesn't require the use of a chest strap, therefore improving ventilation

Secure fixation of the Pad to the OR table with TrenMAX[®] Clamps





Illustrated Complete Set Up: Estape TrenMAX® Pad, 4 each Estape TrenMAX® Clamps

Designed to meet AORN recommendations in "Recommended practices for positioning the patient in the perioperative practice setting" to prevent tissue injury and ischemia that may be caused by tucking a patient's arms at his or her side.

Ordering Information

Description

Cat. No.

Estape TrenMAX® Pad 5/case **415**
 [1- White Base Pad w/8 Rail Straps & The Sticky Pad™, 4- Green Foam Arm Straps]

Estape TrenMAX® Clamp 1 ea. **415-215**
 ~ Estape TrenMAX® Clamps tested and approved ONLY with TrenMAX® Pad ~

Estape TrenMAX® Extensions 20-2 packs/case **415-EXT**
 Biceps Extension Straps for Obese Patients



Designed in conjunction with Ricardo Estape, M.D. Miami FL
 TrenMAX® is a registered trademark of Innovative Medical Products, Inc.
 PATENT APPLIED FOR
 TrenMAX® Clamp - PATENT APPLIED FOR
 The Sticky Pad™ is a trademark of Innovative Medical Products, Inc.
 AORN is a registered trademark of AORN, Inc. AORN does not
 endorse any commercial company's products or services.
 © November 2016 Innovative Medical Products, Inc.
 ALL RIGHTS RESERVED

imp Innovative
 Medical Products, Inc.

87 Spring Lane, PO Box 8028
 Plainville, CT 06062

PH: (860) 793-0391 FAX: (866) 459-1805

Toll Free PH: 1-800-467-4944

www.impmedical.com