


# De Mayo D2<sup>®</sup>


KNEE POSITIONER

## The green standard in Knee Positioning




### KEY FEATURES

- Reduced Procedure Fatigue**  
 The De Mayo D2 Knee Positioner supports the weight of the limb and acts like an extra set of hands, reducing staff fatigue
- Enhanced Stability and Visualization**  
 Creates a more stable operative site with enhanced visualization for all knee replacement surgeries including robotic assisted surgeries
- Easily Adjustable During Surgery**  
 The patented carriage works in conjunction with the baseplate to easily adjust the leg during knee surgery
- Simple Sterile Setup**  
 Sterile baseplate assembly is delivered directly to the OR table
- Easier to Clean, Sterilize, and Store**  
 Knee Positioner can be stored in filtered sterilization containers for convenient processing in SPD



**CAT #923-S**  
(positioner not included)

Accommodates: 25" or 30" De Mayo Knee Positioner base plate for Knee Positioner



**CAT #923-M**  
(positioner not included)

Accommodates: De Mayo Boot, Clamp, and Universal Distractor for Knee Positioner

### ORDERING INFORMATION

Catalog Number	Description
803-114*	De Mayo D2 <sup>®</sup> Knee Positioner <sup>®</sup>
803-GP-10	Sterile Pressure Protector Pad Set
803-ABD	90° Aluminum Boot for De Mayo Knee Positioner
713-717	De Mayo Push Button™ Pin Locking Clamp
120	Quad Clamp <sup>®</sup>
903	Single Lever Clamp
621	De Mayo Universal Distractor <sup>®</sup> 2.0
907-GP-10	Sterile Pressure Protective Pad Set for De Mayo Universal Distractor <sup>®</sup> 2.0

#### Filtered Sterilization Containers for De Mayo D2 Knee Positioner<sup>®</sup>

Catalog Number	Description
923-S	De Mayo Knee Positioner <sup>®</sup> Small Filtered Sterilization Container & Insert (positioner not included)
923-M	De Mayo Knee Positioner <sup>®</sup> Medium Filtered Sterilization Container & Insert (positioner not included)
923-M Kit	De Mayo Knee Positioner <sup>®</sup> Filtered Sterilization Container Kit containing 923-S & 923-M (positioner not included)

\*Various configurations available. Call IMP for a quote and additional information



IMP<sup>®</sup> products are protected by patent & patent pending rights ~ go to IMPpatents.com  
 ©2024 Innovative Medical Products, Inc.  
 ALL RIGHTS RESERVED



Innovative  
 Medical  
 Products™

87 Spring Lane,  
 Plainville, CT 06062  
 PH: (800) 467-4944 F: (866) 459-1805  
 www.IMPmedical.com

SEL-IMP-0007