

NEW

SteriPod™ Filtered Sterilization Containers

Designed to accommodate the De Mayo Knee Positioners® & Accessories

KEY FEATURES

- Eliminate Blue Wrapping**

By using the SteriPod™ Filtered Sterilization Containers, you can eliminate the need for traditional blue wrapping of your positioners. This not only ensures a more efficient process but also alleviates concerns about excess waste and delays caused by potential tears

- Lower Operating Expenses**


Our containers aim to decrease the overall costs associated with the sterilization process by ensuring that medical instruments are properly sterilized and maintained

- Easier Storage**

Stackable containers for easy storage so that you can store multiple containers in a small area, while prolonging the life of your equipment


- Support Green Purchasing**

The use of these containers aligns with environmentally conscious practices in purchasing and disposal by eliminating the use of blue wrap




CAT #923-A2 (positioner not included)
Accommodates: Adapt2Fit (entire system)

Inner Container Dimensions	L 21.675"	W 10.25"	H 9.125"
Outer Container Dimensions	L 23.675"	W 11"	H 9.875"
Weight with Positioner	24.9 lbs.		



CAT #923-S (positioner not included)
Accommodates: 25" or 30" De Mayo Knee Positioner base plate for Knee Positioner

Inner Container Dimensions	L 32.125"	W 9.5"	H 4.4"
Outer Container Dimensions	L 34"	W 10.2"	H 5"
Weight with Positioner	22.85 lbs.		



CAT #923-M (positioner not included)
Accommodates: De Mayo Boot, Clamp, and Universal Distractor for Knee Positioner

Inner Container Dimensions	L 16.8"	W 10.2"	H 7.2"
Outer Container Dimensions	L 18.75"	W 11"	H 7.8"
Weight with Accessories	14.4 lbs.		

*SteriPod™ Filtered Sterilization Containers are available with or without the De Mayo Knee Positioner. Contact IMP for more information.

ORDERING INFORMATION

SteriPod Filtered Sterilization Containers

Catalog Number	Description
923-A2	De Mayo Adapt2Fit® SteriPod™ Filtered Sterilization Container & Insert (positioner not included)
923-A2-CO	De Mayo Adapt2Fit® SteriPod™ Filtered Sterilization Container Only (positioner not included)
923-S	De Mayo Knee Positioner® Short SteriPod™ Filtered Sterilization Container & Insert (positioner not included)
923-M	De Mayo Knee Positioner® Medium SteriPod™ Filtered Sterilization Container & Insert (positioner not included)
923-M KIT	De Mayo Knee Positioner® SteriPod™ Filtered Sterilization Container Kit containing 923-S & 923-M (positioner not included)
CM-10011-851-B	De Mayo Knee Positioner® SteriPod™ System Single Use Paper Filters (100pk)

The filter size is 7 3/8"



IMP® products are protected by patent & patent pending rights - go to IMPpatents.com
©2024 Innovative Medical Products, Inc.
ALL RIGHTS RESERVED



Positioned to Perform

Innovative Medical Products™

87 Spring Lane,
Plainville, CT 06062
PH: (800) 467-4944 F: (866) 459-1805
www.IMPmedical.com

SEL-IMP-0022 V2